

FavorPrepTM
Genomic DNA Mini Kit
(Blood/Cultured Cell)

User Manual

Cat. No.: FABGK 100 (100 Preps)
FABGK 300 (300 Preps)

For Research Use Only

Kit Contents

Cat. No. / preps	FABGK 100 (100 Preps)	FABGK 300 (300 Preps)
RBC Lysis Buffer	135 ml	405 ml
FATG Buffer	30 ml	75 ml
FABG Buffer	40 ml	100 ml
W1 Buffer	45 ml	130 ml
Wash Buffer** (concentrated)	25 ml	50 ml
Elution Buffer	30 ml	75 ml
FABG Column	100 pcs	300 pcs
2 ml Collection Tube	200 pcs	600 pcs

**Add 100 ml / 200 ml of ethanol (96~100%) to Wash Buffer when first open.

Fresh Blood Protocol:

• Step 1 – RBC Lysis

1. Collect fresh human blood in an anticoagulant-treated collection tube.
2. Transfer up to 300µl fresh blood to a microcentrifuge tube (not provided).
If the blood sample is more than 300 µl (up to 1 ml), add the sample to a sterile 15 ml centrifuge tube.
3. Add 3 x the sample volume of RBC Lysis Buffer and mix by inversion. Do not vortex.
4. Incubate at room temperature for 10 minutes.
5. Centrifuge at 3000 x g for 5 minutes and completely remove the supernatant.
6. Resuspend the pellet with 100 µl of RBC Lysis Buffer.

Step 2 – Cell Lysis

7. Add 200 µl of FABG Buffer and mix by vortex.
8. Incubate for 10 minutes at room temperature or until the sample lysate is clear. During incubation, invert the tube every 3 minutes.
9. Briefly spin the tube to remove drops from the inside of the lid.
10. Preheat Elution Buffer (200 µl per sample) in a 70°C water bath (for Step 5 DNA Elution).
11. (Optional Step): If RNA-free genomic DNA is required, add 5µl of 10 mg/ml RNase A to the sample and vortex, then incubate for 5 minutes at room temperature.

Step 3 – Binding

12. Add 200µl ethanol (96~100%) to the sample. Mix thoroughly by pulse-vortexing for 10 seconds. (Pipetting if there is any precipitate.)
13. Briefly spin the tube to remove drops from the inside of the lid.
14. Place a FABG Column to a collection tube. Transfer the sample mixture (including any precipitate) carefully to FABG Column. Centrifuge for 5 minute and discard the flow-through then place FABG Column to a new Collection tube.

Step 4 – Washing

15. Immediately, Wash FABG Column with 400µl W1 Buffer (ethanol added) by centrifuge for 30 seconds then discard the flow-through.
 - Make sure that ethanol has been added into W1 Buffer when first open.
16. Wash FABG Column with 600µl Wash Buffer (ethanol added) by centrifuge for 30 seconds then discard the flow-through.
 - Make sure that ethanol has been added into Wash Buffer when first open.
17. Centrifuge for an additional 3 min to dry the column.
 - Important Step! This step will avoid the residual liquid to inhibit subsequent enzymatic reactions.

Fresh Blood Protocol

Step 5 – Elution

18. Place FABG Column to a new 1.5ml microcentrifuge tube.
19. Add 100µl of Preheated Elution Buffer or TE to the membrane center of FABG Column. Stand FABG Column for 3~5 min or until the buffer is absorbed by the membrane.
 - Important Step! For effective elution, make sure that the elution solution is dispensed onto the membrane center and absorbed completely.
20. Centrifuge for 30 seconds to elute the DNA .
 - Standard volume for elution is 100 µl. If sample has low number of cells, reduce the elution volume (30 µl - 50 µl) to increase DNA concentration. If higher DNA yield is required, repeat the DNA Elution step to increase DNA recovery and the total volume could be 200 µl.

Step Final - Pure DNA

21. Store the DNA fragment at 4°C or -20°C.

Frozen Blood Protocol:

Sample Preparation

1. Transfer up to 200µl sample to a microcentrifuge tube (not provided). If the sample volume is less than 200µl, add the appropriate volume of PBS.
2. Add 30µl Proteinase K (10 mg/ml) to the sample. Mix thoroughly by pulse-vortexing. Incubate for 15 minutes at 60°C.

Step 1 –Cell Lysis

3. Add 200µl FABG Buffer to the sample. Mix thoroughly by pulse-vortexing.
4. Incubate at 70°C for 15 minutes to lyse the sample. During incubation, invert the sample every 3~5 minutes.
5. Briefly spin the tube to remove drops from the inside of the lid.
6. Preheat Elution Buffer (200 µl per sample) in a 70°C water bath (for Step 4 DNA Elution).
7. (Optional Step): If RNA-free genomic DNA is required, add 5µl of 10 mg/ml RNase A to the sample and vortex, then incubate for 5 minutes at room temperature.
8. Follow the Fresh Blood Protocol starting from Step 3 (Binding).

Buffy Coat Protocol:

Sample Preparation

Centrifuge whole blood at 3,300xg for 10 minutes at room temperature and you will get three different fractions: the upper clear layer is plasma; the intermediate layer is buffy coat, containing concentrated leukocytes; the bottom layer contains concentrated erythrocytes. Extraction total DNA from buffy coat will yield 5-10 times more DNA than an equivalent volume of whole blood.

Step 1 – RBC Lysis

1. Transfer up to 200µl buffy coat to a microcentrifuge tube (not provided).
2. Add 3 x the sample volume of RBC Lysis Buffer and mix by inversion.
3. Incubate at room temperature for 10 minutes. During incubation, invert the tube every 3 minutes.
4. Centrifuge for 1 minutes at full speed (14,000 rpm or 10,000 x g) and completely remove the supernatant.
5. Resuspend the pellet with 500 µl of RBC Lysis Buffer, then centrifuge for 1 minutes and completely remove the supernatant.
6. Resuspend the pellet with 200 µl of RBC Lysis Buffer. (Mix the tube by vortex only. Be sure the pellet is completely resuspended or the column would be barred when binding.)

Step 2 – Cell Lysis

7. Add 250 µl of FABG Buffer and mix by vortex.
8. Incubate for 30 minutes at room temperature or until the sample lysate is clear. During incubation, invert the tube every 3 minutes.
9. Briefly spin the tube to remove drops from the inside of the lid.
10. Preheat Elution Buffer (200 µl per sample) in a 70°C water bath (for Step 5 DNA Elution).
11. (Optional Step): If RNA-free genomic DNA is required, add 5µl of 10 mg/ml RNase A to the sample and vortex, then incubate for 5 minutes at room temperature.

Step 3 – Binding

12. Add 250µl ethanol (96~100%) to the sample. Mix thoroughly by pulse-vortexing for 10 seconds. (Pipetting if there is any precipitate.)
13. Briefly spin the tube to remove drops from the inside of the lid.
14. Place a FABG Column to a collection tube. Transfer the sample mixture (including any precipitate) carefully to FABG Column. Centrifuge for 5 minute and discard the flow-through then place FABG Column to a new Collection tube.

Buffy Coat Protocol:

Step 4 – Washing

15. Immediately, Wash FABG Column with 400µl W1 Buffer (ethanol added) by centrifuge for 1 minute then discard the flow-through.
 - Make sure that ethanol has been added into W1 Buffer when first open.
16. Wash FABG Column with 600µl Wash Buffer (ethanol added) by centrifuge for 1 minute then discard the flow-through.
 - Make sure that ethanol has been added into Wash Buffer when first open.
17. Centrifuge for an additional 3 min to dry the column.
 - Important Step! This step will avoid the residual liquid to inhibit subsequent enzymatic reactions.

Step 5 – Elution

18. Place FABG Column to a new 1.5ml microcentrifuge tube.
19. Add 100µl of Preheated Elution Buffer or TE to the membrane center of FABG Column. Stand FABG Column for 3~5 min or until the buffer is absorbed by the membrane.
 - Important Step! For effective elution, make sure that the elution solution is dispensed onto the membrane center and absorbed completely.
20. Incubate the FABG Column at 37°C for 10 minutes in an incubator.
21. Centrifuge for 1 minute to elute the DNA .
 - Standard volume for elution is 100 µl. If sample has low number of cells, reduce the elution volume (30 µl - 50 µl) to increase DNA concentration. If higher DNA yield is required, repeat the DNA Elution step to increase DNA recovery and the total volume could be 200 µl.

Step Final - Pure DNA

22. Store the DNA fragment at 4°C or -20°C.

Cultured Cell Protocol:

Sample Preparation

For Cultured Cells:

- i. Trypsinize the adherent cells before harvesting.
- ii. Transfer the appropriate number of cell (up to 1×10^7) to a 1.5ml microcentrifuge tube (not provided) and centrifuge at $6000 \times g$ for 20 seconds.
- iii. Remove the supernatant and resuspend the cells with 150 μ l of RBC Lysis Buffer.

For Fresh blood (except human blood):

- i. The sample volume of mammalian blood (non-nucleated) can be up to 50 μ l; the sample volume of nucleated erythrocytes (eg. bird or fish) can be up to 10 μ l.
- ii. Add 150 μ l of FATG Buffer and the blood sample into a 1.5ml microcentrifuge tube (not provided), then mix by vortex.

Step 1 –Cell Lysis

1. Add 200 μ l of FABG Buffer and mix by vortex for 5 seconds.
2. Incubate for 10 minutes at 70°C or until the sample lysate is clear. During incubation, invert the tube every 3 minutes.
3. Briefly spin the tube to remove drops from the inside of the lid.
4. Preheat Elution Buffer (200 μ l per sample) in a 70°C water bath. (for Step 4 DNA Elution).
5. (Optional Step): If RNA-free genomic DNA is required, add 5 μ l of 10 mg/ml RNase A to the sample and vortex, then incubate for 5 minutes at room temperature.
6. Follow the Fresh Blood Protocol starting from Step 3 (Binding).

Cultured Cell Protocol:

Sample Preparation

- i. Harvest appropriate number of cell (up to 5×10^7) to a 1.5ml microcentrifuge tube (not provided) and centrifuge at $5000 \times g$ for 10 minute, then discard the supernatant.
- ii. Add 600 μ l of sorbitol buffer (1.2 M sorbitol; 10 mM CaCl₂; 0.1 M Tris-HCl pH 7.5; 35mM mercaptoethanol) and resuspend the pellet.
- iii. Add 200 U of lyticase or zymolase. Incubate for 30 minutes at 30°C.
- iv. Centrifuge the mixture at $2,000 \times g$ for 10 minutes to harvest the spheroplast.
- v. Remove the supernatant and add 200 μ l of FATG Buffer to the tube and resuspend the cell pellet by vortex.
- vi. Incubate at room temperature for 5 minutes.
- vii. Follow the Cultured Cell Protocol starting from Step 1 (Cell Lysis).