

**FavorPrep<sup>TM</sup>**  
**Soil DNA Isolation Mini Kit**

**User Manual**

**Cat. No.: FASOI 001 (50 Preps)**  
**FASOI 001-1 (100 Preps)**

**For Research Use Only**

# Introduction

The FavorPrep™ Soil DNA Isolation Kit is designed for isolation of high quality total DNA from soil sample. The inhibitors of the downstream application such as polysaccharides, humic acid, phenolic compounds will be removed by utilizing the DNA binding column and the buffer system. The entire procedure is not required the phenol-chloroform extraction and can be finished within 60 minutes. The purified DNA is ready for PCR and other downstream application.

# Specification:

Sample Size: 0.2-1 g of soil sample
Handling Time: about 60 minutes

# Kit Contents

	FASOI 001 (50 preps)	FASOI 001-1 (100 preps)
Glass Beads	12 g	25 g
SDE1 Buffer	40 ml	70 ml
SDE2 Buffer	15 ml	25 ml
SDE3 Buffer	15 ml	30 ml
SDE4 Buffer	25 ml	40 ml
Wash Buffer	20 ml*	40 ml**
Elution Buffer	32 ml	35 ml
SDE Mini Column	50 pcs	100 pcs
1.5 ml Tube (for Elution)	50 pcs	100 pcs
2.0 ml Tube (for Collection)	100 pcs	200 pcs
Bead Tube	50 pcs	100 pcs


\*Add 80 ml of ethanol (96~100%) to Wash Buffer when first open.

\*\*Add 140 ml of ethanol (96~100%) to Wash Buffer when first open.

## Important Notes:

1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
2. Check SDE1 Buffer before use, Warm SDE1 Buffer at 60°C for 10 minutes if any precipitate formed.
3. Add 80 ml/ 140 ml of ethanol (96-100%) to Wash Buffer when first open.
4. Prepare a water bath to 70 °C before the operation.
5. All centrifuge steps are done at full speed (14,000 rpm or 10,000 x g) in a microcentrifuge.
6. Preheat Elution Buffer or ddH<sub>2</sub>O to 60°C for elution step.

## Brief Procedure:



# General Protocol:

## Please Read Important Notes Before Starting Following Steps.

- 1. Add 200 mg of Glass Beads into a 2.0 ml Bead Tube (provided). And transfer 0.25~1 g of soil sample into Bead Tube then place on ice.**
  - If the sample is liquid, add 200  $\mu$ l of sample into a 2.0 ml Beads Tube.
- 2. Add 600  $\mu$ l of SDE1 Buffer to the sample, vortex at maximum speed for 5 minutes. Incubate the sample at 70 °C for 10 minutes and vortex the sample twice during the incubation.**
  - For isolation of DNA from gram positive bacteria, do a further incubation at 95 °C for 5 minutes.
- 3. Briefly spin the tube to remove drops from the inside of the lid.**
- 4. Cool down the sample and add 200  $\mu$ l of SDE2 Buffer to the sample, mix well by vortexing. Incubate the sample on ice for 5 minutes.**
- 5. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 5 minutes.**
- 6. Carefully transfer the clarified lysate to a 1.5 ml microcentrifuge tube (not provided). And measure the volume of the clarified lysate.**
  - Avoid pipetting any debris and pellet.
- 7. Add 1 volume of isopropanol, vortex to mix well. centrifuge at full speed (14,000 rpm or 10,000 x g) for 10 min to pellet DNA.**
  - For example: If the clarified lysate volume is 450  $\mu$ l, add 450  $\mu$ l of isopropanol to the clarified lysate.
- 8. Carefully discard the supernatant and invert the tube on the paper towel for 1 min to remove residual liquid.**
  - Do not disturb the pellet.
- 9. Add 200  $\mu$ l of pre-heated Elution Buffer or ddH<sub>2</sub>O, vortex to dissolve the DNA pellet completely.**
- 10. Add 100  $\mu$ l of SDE3 Buffer to the sample, mix well by vortexing. Incubate the sample at room temperature for 2 minutes.**
  - Note:** SDE3 Buffer must be suspended completely by vigorously vortexing before every using.
  - use 1ml pipettor and cut off the end of 1 ml tip to make it easier for pipetting the SDE3 Buffer.

11. **Centrifuge at full speed (14,000 rpm or 10,000 x g) for 2 minutes.**
12. **Carefully transfer the clarified lysate to a 1.5 ml microcentrifuge (not provided). And measure the volume of the clarified lysate.**  
--Avoid pipetting any debris and pellet.
13. **(Optional) If RNA-free DNA is required, add 1  $\mu$ l of 100 mg/ml RNase A (not provided) to the sample and mix well. Incubate at room temperature for 2 min.**
14. **Briefly spin the tube to remove drops from the inside of the lid.**
15. **Add 1 volume of SDE4 Buffer and 1 volume of ethanol (96~100%) to the clarified lysate, mix thoroughly by pulse-vortexing.**  
For example: If the clarified lysate volume is 250  $\mu$ l, add 250  $\mu$ l of SDE4 Buffer and 250  $\mu$ l of ethanol (96~100%) to the clarified lysate.
16. **Place a SDE Column into a Collection Tube and transfer all of the sample mixture to the SDE Column. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min then discard the flow-through. Place the SDE Column to a new Collection Tube.**
17. **Add 750  $\mu$ l of Wash Buffer (ethanol added) to SDE Column. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min then discard the flow-through. And repeat this step for one more time.**  
--Make sure that ethanol (96~100%) has been added into Wash Buffer when first open.
18. **Centrifuge at full speed (14,000 rpm) for an additional 3 min to dry the SDE column.**  
--**Important step!** This step will avoid the residual liquid to inhibit subsequent enzymatic reactions.
19. **Place SDE Column into a Elution Tube, Add 50~200  $\mu$ l of preheated Elution Buffer or ddH<sub>2</sub>O to the membrane center of the SDE Column. Stand the SDE Column for 2 min at room temperature.**  
--**Important step!** For effective elution, make sure that the Elution Buffer or ddH<sub>2</sub>O is dispensed onto the membrane center and is absorbed completely.
20. **Centrifuge for 1 min to elute DNA.**

## Troubleshooting:

### Low yield

- Too many sample were used
  - reduce the sample volume.
- Poor sample lysis because of insufficient sample beating with glass beads
  - Extend the vortexing time at maximum speed
- Poor sample lysis because of insufficient incubation time
  - Extend the incubation time
- Ethanol is not added into the lysate before transferring into SDE Mini Column
- Ethanol is not added into Wash Buffer when first open; the volume or the percentage of ethanol is not correct before adding into Wash Buffer
- Elution of DNA is not efficient
  - Make sure the pH of ddH<sub>2</sub>O is between 7.5- 8.5.
  - After Elution Buffer or ddH<sub>2</sub>O is added, stand the SDE Mini Column for 2 min before centrifugation.

### Column is clogged

- Disrupt the debris pellet when transfer the sample supernatant
- Sample is too viscous
  - Reduce the sample volume.

### Purified DNA dose not perform well in downstream application

- Sample is old
  - Always use fresh or well-stored sample for stool DNA extraction.
- Residual ethanol contamination
  - After Wash step, centrifuge at 4,000 x g for an additional 3 minutes to dry the SDE Mini Column.
- RNA contamination