

# FavorPrep<sup>TM</sup> Plant Total RNA Mini Kit (for woody plant)


Cat.: FAPRK 003 (50Preps)  
FAPRK 003-1 (100Preps)  
(For Research Use Only)

## Kit Contents

	FAPRK 003	FAPRK 003-1
FARB-1 Buffer	30 ml	60 ml
FARB-2 Buffer	4 ml	8 ml
Wash Buffer 1	30 ml	60 ml
Wash Buffer 2 (conc.)*	15 ml	35 ml
RNase-free water	6 ml	6 ml
FARB Mini Column	50 Pcs	100 Pcs
2 ml Collection Tube	50 Pcs	100 Pcs

\*Add 60 / 140 ml ethanol (96~100%) to Wash Buffer when first open.

## Brief Procedure



## Important Notes

1. Make sure everything is RNase-free when handling RNA.
2. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
3. Pipet a required volume of FARB-1 Buffer to another RNase-free container and add 10 ul of B-mercaptoethanol (B-ME) per 1ml FARB-1 Buffer before use.
4. Add required amount of ethanol(96-100%) as bottle indicated to Wash Buffer 2 when first open.
5. Dilute RNase-free DNase 1 in reaction buffer (1M NaCl, 10mM MnCl<sub>2</sub>, 20mM Tris-HCl, pH7.0 at 25 °C) to final conc.= 0.5U/ul.

## General Protocol:

Please Read Important Notes Before Starting Following Steps.

1. Grind up to 100 mg plant sample under liquid nitrogen to a fine powder and transfer to a new microcentrifuge tube (not provided).
2. Add 500 ul of FARB-1 Buffer (B-ME added) to the sample powder and vortex vigorously.
3. Add 50 ul of FARB-2 Buffer and incubate at 70 °C for 10 min, vortex every 3 min during incubation.
4. Centrifuge at 12,000 rpm for 5 min at 5 °C.
5. Transfer the clarified supernatant to a new microcentrifuge tube (not provided) and adjust the volume of the clear lysate.  
--Avoid pipett any debris and pellet from the Collection Tube.

**6. Add 0.9 volume of ethanol (96-100%) to the clear lysate and mix well.**

--For example, add 450 ul of ethanol to 500 ul of clear lysate.

**7. Place a FARB Mini Column into a Collection Tube, And transfer 750 ul of the ethanol added sample mixture (including any precipitate) to FARB Mini Column. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min and discard the flow-through.**

**8. Repeat step 7 for rest of the sample mixture.**

**9. (Optional): To eliminate genomic DNA contamination of RNA, follow the steps from 9a. Otherwise, proceed to step 10 directly.**

**9a. Add 250 ul of Wash Buffer 1 to FARB Mini Column, Centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min then discard the flow-through.**

**9b. Add 100 ul of RNase-free DNase 1 solution (0.5 U/ul, not provided) to the membrane center of FARB Mini Column. Place the Column on the benchtop for 15 min.**

**9c. Add 250 ul of Wash Buffer 1 to FARB Mini Column. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min then discard the flow-through.**

**9d. After DNase 1 treatment, proceed to step 11.**

**10. Add 500 ul of Wash Buffer 1 to wash FARB Mini Column, Centrifuge for 1 min then discard the flow-through.**

**11. Wash FARB Mini Column twice with 750 ul of Wash Buffer 2 by centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min then discard the flow-through.**

--Make sure that ethanol has been added into Wash Buffer 2 when first open.

**12. Centrifuge at full speed (14,000 rpm or 10,000 x g) for an additional 3 min to dry the column.**

--Important Step! This step will avoid the residual liquid to inhibit subsequent enzymatic reaction.

**13. Place FARB Mini Column to Elution Tube.**

**14. Add 50 ul of RNase-free ddH<sub>2</sub>O to the membrane center of FARB Mini Column. Stand FARB Mini Column for 1 min.**

--Important Step! For effective elution, make sure that the elution solution is dispensed of the membrane center and is absorbed completely.

**15. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 2 min to elute RNA.**

**16. Store RNA at -70 °C.**