SuperPicTure™ Polymer

Detection Kit

Zymed® Polymer Detection System

ZYMED® Laboratories

invitrogen immunodetection

<u>CAT. NO.</u>	<u>FOR PRIMARY ANTIBODY</u>	<u>ENZYME</u>	<u>GOOD FOR</u>
87-8963	Broad Spectrum*	HRP Polymer**	1000 slides

*Will react with mouse, rabbit, guinea pig and rat primary antibodies.

**This kit does not contain a chromogen. AEC *Single Solution* (Cat. No. 00-1111) or DAB Kit (Cat. No. 00-2014) are recommended to be used with this kit.

INTENDED USE

For In Vitro Diagnostic Use. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

SuperPicTure™ Polymer Detection kits are designed to reveal antigens that have reacted with a user-supplied primary antibody on human tissue or cell samples. Frozen or paraffin-embedded tissues, freshly prepared lymphocytes, and fixed cultured cells can be used.

BACKGROUND

The *SuperPicTure™* Polymer Detection kit features Zymed's proprietary and enhanced HRP polymer technology to provide excellent sensitivity, high specificity, and a shorter protocol than conventional immunohistochemical staining methods. Zymed's new polymer conjugation methods have resulted in a versatile polymer capable of intense nuclear, cytoplasmic, and membrane antigen staining. Drawbacks of other polymer systems include compromised cell penetration ability that can cause weak nuclear staining.¹ Zymed's polymer conjugate contains multiple HRPs to provide excellent sensitivity and specificity. The polymer staining method is faster and easier than a conventional LAB-SA protocol because it uses an HRP polymer conjugated to the second antibody and eliminates one incubation step. The *SuperPicTure™* Polymer detection system is also cleaner because it uses only one protein-based reagent, whereas the conventional LAB-SA method uses two protein-based reagents. In addition, *SuperPicTure™* Polymer kits do not contain biotin nor streptavidin, so potential background due to endogenous biotin activity is completely avoided. Another advantage of Zymed's *SuperPicTure™* Polymer detection system is the shorter polymer incubation time; the *SuperPicTure™* system only requires 10 minutes, while leading competitors' polymers require 30 minutes.

Zymed offers prediluted *2ndGen* primary antibodies that are ready-to-use with *SuperPicTure™* Polymer Detection kits. For a complete list of our primary antibodies, call for a free copy of Zymed Pathology Catalog or contact Zymed's technical service for details.

PROTOCOL SUMMARY

After formalin-fixed, paraffin-embedded tissue sections are deparaffinized in xylene and dehydrated in a graded series of ethanol, an endogenous peroxidase quenching step is performed. This endogenous peroxidase activity can be quenched by either using 3% hydrogen peroxide in methanol or Zymed's Peroxo-Block™ (Cat. No. 00-2015). No blocking solution is required. After incubating the primary antibody on the tissue, the HRP polymer is added, and after a wash step, the chromogen is then added, and the peroxidase catalyzes the substrate (hydrogen peroxide) and converts the chromogen to a colored deposit, which visualizes the location of the antigen.

STORAGE AND SHELF LIFE

Store kit at 2-8°C. All performance claims are void after kit expiration date.

REAGENTS SUPPLIED

Reagent A. One dropper bottle (110 mL) of ready-to-use HRP polymer conjugate.

REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

- Xylene, ethanol, and absolute methanol
- Distilled or deionized water
- 30% hydrogen peroxide
- 10 mM phosphate-buffered saline (PBS) (Optional: 0.05% Tween 20)
- Mini PAP Pen (Zymed Cat. No 00-8877)
- Primary antibody
- HRP Substrate/Chromogen reagents:
 - AEC *Single Solution* (Zymed Cat. No. 00-1111)
 - DAB Kit (Zymed Cat. No. 00-2014)
 - DAB-Black™ Kit (Zymed Cat. No. 00-2013)
- Hematoxylin (Zymed Cat. No. 00-8001)
- Mounting Media
 - Clearmount™ (Zymed Cat. No. 00-8110) or
 - Histomount™ (Zymed Cat. No.00-8030)

www.invitrogen.com

Invitrogen Corporation • 1600 Faraday Avenue • Carlsbad • CA 92008 • Tel: 800.955.6288 • E-mail: tech_support@invitrogen.com

PIN: 30756

Rev. 12/06

DCC-06-0084

Copyright 2005, Zymed Laboratories

SUGGESTED STAINING PROTOCOL

A. PRELIMINARY PREPARATION OF SLIDES

<u>1. SPECIMEN PREPARATION</u>	<p>-Appropriate tissue and antigen fixation is required to obtain reproducible performance and reliable interpretations. Suitable fixatives include: 10% neutral buffered formalin, B5, Bouin's, Zinc formalin or alcohol-base fixatives.</p> <p>-10% neutral buffered formalin, B5, Bouin's, Zinc formalin or alcohol-base fixatives are regarded as suitable fixatives for most antigens of clinical significance.</p> <p>- Formalin-fixed tissues post-fixed in B5 before paraffin embedding may show improved stain. Cell smears prepared from body fluids should be made to assure a monolayer of cells. Multilayers of cells can trap staining reagents and interfere with the interpretation of the results. Smears should be fixed immediately after preparation. Depending on the properties of the antigen, cell smears are usually stable for one to two weeks when stored at 4°C.</p> <p>- Fixation of cytospin or frozen sections can be done with acetone (100%) at 4°C for a period of 10 minutes.</p>
<u>2. SLIDE PREPARATION</u>	<p>-Precoat slides with HistoGrip™ (Zymed Cat. No. 00-8050). As an alternative, precoat with 0.1% poly-L-lysine in water, then air dry.</p>
<u>3. DEPARAFFINIZATION AND REHYDRATION</u>	<p>-Paraffin sections are deparaffinized with xylene and rehydrated in a graded series of ethanol. Wash cell smears or tissue in a PBS bath for 10 minutes before starting the staining procedure. Note: Do not allow specimen or tissue sections to dry from this point on.</p>
<u>4. CONTROL SLIDES</u>	<p>Three control slides are necessary for the interpretation of results.</p> <ul style="list-style-type: none">• Positive Tissue Control: A specimen, processed in the same way as the unknown, contains the antigen to be stained.• Reagent Control/Negative Control-1: An additional slide that will be treated with a non-immune serum instead of same concentration of primary antibody. Any staining observed on the specimen is probably due to non-specific protein binding or non-specific binding of other reagents.• Negative Control-2: A specimen, processed in the same way as the unknown, does not contain the antigen to be stained [optional].

SUGGESTED STAINING PROTOCOL CONTINUED

B. STAINING PROTOCOL

Reagent	Staining Procedures	Incubation Time (Min.)
<u>1. PEROXIDASE QUENCHING SOLUTION</u> (Not provided)	<p>This step is optional. Do step 1 if elimination of endogenous peroxidase activity is necessary.</p> <p>-For paraffin-embedded tissues,</p> <ol style="list-style-type: none">Add 1 part of 30% hydrogen peroxide to 9 parts of absolute methanol. Mix well.Submerge slides in PEROXIDASE QUENCHING SOLUTION for 10 minutes.Wash with PBS 2 min., 3 times. Proceed to Step 2. <p>For frozen tissue, treat with Peroxo-Block™ (cat. no. 00-2015) for 45 sec. Wash immediately. Peroxo-Block™ can also be used with paraffin tissues.</p>	10 (optional)
<u>2. PRIMARY ANTIBODY</u> (Not provided)	<p>Optimal dilution and incubation times should be determined by the investigator. 2nd Gen Primary antibodies are optimized for use with SuperPicTure™ detection kits. (Dilution and incubation time are dependent on sample preparation, antibody affinity, amt. of antigen present, and antigen accessibility).</p> <ol style="list-style-type: none">Apply 2 drops (100 µL) or enough to completely cover tissue, of PRIMARY ANTIBODY to each section.Incubate in moist chamber for 30 min.Rinse with PBS for 2 min., 3 times.	30-60
<u>3. HRP POLYMER CONJUGATE</u> Ready-to-use. Reagent A	<ol style="list-style-type: none">Apply 2 drops (100 µL) or enough to completely cover tissue, of HRP POLYMER CONJUGATE to each section. Incubate for 10 min.Rinse with PBS for 2 min., 3 times.	10
<u>4. CHROMOGEN</u> (Not included)	<ol style="list-style-type: none">Prepare and add chromogen by following the manufacturer's instructions. We recommend using either the AEC Single Solution (Cat. No. 00-1111) or DAB Kit (Cat. No. 00-2014).	DAB 5 AEC 10

From this point on, the investigator should proceed with his or her own established laboratory protocols for counterstaining and mounting. Zymed recommends using Hematoxylin (Cat. No. 00-8001) as a counterstain and either Histomount™ (Cat. No. 00-8030) or Clearmount™ (Cat. No. 00-8010) as a mounting solution for DAB.

INTERPRETATION OF RESULTS

The following table interprets results which might be achieved using various staining controls.

Case No.	+ control	reagent control	(-) control	Sample	Analysis
1	-	-	-	-	Procedure performed incorrectly.
2	+	+	+	+	Non-specific staining due to protein binding or endogenous peroxidase activity.
3	+	-	+	+/-	Negative control contains the antigen.
4	-	-	-	+	Positive control does not contain the antigen.
5	+	-	-	-	Sample does not contain the antigen.
6	+	-	-	+	Sample contains the antigen.

TROUBLESHOOTING

Possible causes for negative staining on positive slides:

1. Steps in the staining protocol were performed in incorrect sequence.
2. Primary or secondary antibody incubation steps were omitted.
3. Labile antigens were destroyed.
4. Specimen was improperly fixed and/or processed.
5. Specimen dehydrated during staining.

Possible causes for weak staining on all slides:

1. Specimen retained excess liquid after rinsing steps.
2. Incubation times were insufficient.
3. Substrate prepared improperly.
4. Deparaffinization was incomplete (staining may be accompanied by high background).

Possible causes for high background staining:

1. Endogenous peroxidase activity was incompletely blocked.
2. Deparaffinization was incomplete.
3. Excessive application of tissue adhesive.
4. Inadequate rinsing of slides.
5. Over-development of substrate.
6. Dehydration of specimen during staining.

REFERENCES

1. Petrosyan K, et al. *J Histotechnology* 25(4): 247-250, 2002.

TRADEMARKS

SuperPicTure™, *CAS-Block™*, *Clearmount™*, *HistoGrip™*, *Histomount™*, and *Peroxo-Block™* are trademarks of Zymed Laboratories, Inc. Zymed® and *Histostain®* are registered trademarks of Zymed Laboratories, Inc.

Authorized Representative for IVDD 98/79/EC: Invitrogen Ltd.
Inchinnan Business Park
3 Fountain Drive
Paisley
PA4 9RF
UK

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